



It's a simple
MUST...

**Fit, Function And
Esthetics**

“Keeping your patients smiling”

ClassiTemps

THE NEW STANDARD

ClassiTemps are the new standard in provisional temporization. With the high level of patient education and demands for esthetic dentistry, it's important that in the provisional stages of restorative treatment, the esthetic and functional requirements are satisfied.

STRENGTH AND DURABILITY

ClassiTemps are strong and durable by nature due to the heat and pressure processing. For longer span bridges a wire or cast reinforcement is added.

SIMPLICITY

After your preparations, you simply relin and cement the **ClassiTemps**.



Mitech Dental Laboratory Inc.

5230 South Service Road
Burlington, ON L7L 5K2

Tel: (905) 333-0765

Fax: (905) 333-3763

Toll Free: (800) 509-8815

www.mitechdental.com

info@mitechdental.com

*For decreased chair times, ask about seating
stents. If you would like more detailed
information, please contact us today.*